

Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

Safety Data Sheet

Section 1: Identification				
1.1	Product Identifiers Product Name	1000ppm Ammonium Standard Calibration Solution		
1.2	Application Intended Use	Calibration of ammonium ISE		
1.3	Details of Supplier Manufacturer/Supplier	Sentek Ltd Unit 6-7 Crittall Court Crittall Drive Springwood Industrial Estate Braintree Essex CM7 2SE England		
	Telephone Fax E-mail	+44 (0) 1376 340456 +44 (0) 1376 340453 <u>sentekuk@btconnect.com</u>		
1.4	Emergency Telephone			

Emergency Contact +44

+44 (0) 1376 340456

Section 2: Hazards Identification

Classification of the Substance or Mixture
Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008



Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate, Braintree, Essex, CM7 2SE

Registered in England No. 2753857

Registered Office:

BS EN ISO 9001:2015 Certificate No. FM 36099



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

2.2 Label Elements None. Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008

2.3 Other Hazards None

Section 3: Composition/Information on Ingredients

3.1 Substances

Name	Ammonium Chloride
Formula	NH4Cl
CAS No	12125-02-9
Concentration	<0.5%

Name	Water
Formula	H ₂ O
CAS No	7732-18-5
Concentration	>99.5%

Section 4: First Aid Measures

4.1 Description of First Aid Measures

In Case of Eye Contact	Flush eyes thoroughly with water
In Case of Skin Contact	Rinse skin with water
In Case of Swallowing	Rinse mouth with water. Consult a physician

Section 5: Firefighting Measures

5.1 Extinguishing Media

Suitable Extinguishing Media Water Spray, Foam, Dry Powder, Carbon Dioxide



BS EN ISO 9001:2015

Certificate No. FM 36099

Registered Office:



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

- 5.2 Special Hazards Arising From the Substance or Mixture Hazardous Combustion Products At high temperatures potential release of hydrogen chloride gas, nitrogen oxides and ammonia
- 5.3 Advice for firefighters Wear breathing apparatus if necessary

Section 6: Accidental Release Measures

- 6.1 Personal Precautions, Protective Equipment and Emergency Procedures Use suitable laboratory precautions
- 6.2 Environmental Precautions Do not allow to enter drains
- 6.3 Methods and Materials for Containment and Cleaning up Use suitable absorbent material and dispose of safely

Section 7: Handling and Storage

- 7.1 Precautions for Safe Handling No special measures required other than normal laboratory precautions
- **7.2 Conditions for Safe Storage, Including any Incompatibilities** Keep container tightly closed in a dry and well-ventilated place.
- **7.3** Specific End Uses No other uses other than those specified in Section 1.2

Section 8: Exposure Controls/Personal Protection

8.1 Control Parameters

Components with Workplace Control Parameters None



BS EN ISO 9001:2015 Certificate No. FM 36099

Registered Office: Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate, Braintree, Essex, CM7 2SE



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

8.2 Exposure Controls

General Protective and Hygienic Measures

Usual precautionary laboratory methods for handling chemicals should be followed

Personal Protective Equipment

Eye Protection Use suitable laboratory standard eye protection Skin Protection Handle with gloves Body Protection Not required Respiratory Protection Not required

Environmental Exposure Controls

Do not allow to enter drains

Section 9: Physical and Chemical Properties

9.1 Information on Physical and Chemical Properties

Form:	Liquid
Colour:	Colourless
Odour:	Odourless
Odour Threshold:	Not determined
рН:	Not determined
Melting Point/	Not determined
Freezing Point:	
Boiling Point/	100°C
Boiling Range:	
Flash Point:	Not applicable
Flammability:	Not applicable
Auto-ignition	Not applicable
Temperature:	



BS EN ISO 9001:2015 Certificate No. FM 36099

Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate, Braintree, Essex, CM7 2SE

Registered in England No. 2753857

Registered Office:



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

Decomposition Temperature: Explosion limits: Vapour Pressure: Density: Vapour Density: Evaporation Rate: Water Solubility: Partition Coefficient: n-octanol/water Viscosity: Explosive Properties: Oxidizing Properties:

Not applicable Not determined Not determined Not determined Not determined Soluble Not determined

Not determined

Not determined Not determined Not determined

9.2 Other Information

No further relevant information available

Section 10: Stability and Reactivity

- **10.1 Reactivity** No relevant information available
- **10.2 Chemical Stability** Stable under normal conditions
- **10.3 Possibility of Hazardous Reactions** None known
- **10.4 Conditions to Avoid** No relevant information available
- **10.5** Incompatible Materials Strong oxidizing reagents, strong acids, strong bases

10.6 Hazardous Decomposition products In the event of fire: see Section 5



BS EN ISO 9001:2015 Certificate No. FM 36099

Registered Office: Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate, Braintree, Essex, CM7 2SE



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

Section 11: Toxicological information

- 11.1 Information on Toxicological Effects
 - Ammonium Chloride: Category 4 (H302) Acute Toxicity: **Skin Irritation:** No data available Eye damage/Irritation Eye irritation Sensitization No data available Germ Cell Mutagenicity: No data available **Carcinogenicity:** No data available **Reproductive Toxicity:** No data available Specific Target Organ No data available Toxicity No data available **Aspiration Hazard:**

Section 12: Ecological information

- **12.1 Toxicity** Ammonium chloride – Cyprinus carpio (Carp) (LC50: 209mg/L/96h) Water fleas (EC50: 101mg/L/48h)
- 12.2 Persistence and Degradability No data available
- 12.3 Bio-accumulative Potential No data available
- 12.4 Mobility in Soil No data available
- **12.5** Result of PBT and vPvB No data available
- 12.6 Other Adverse Effects No data available



BS EN ISO 9001:2015 Certificate No. FM 36099

Registered Office:



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

Section 13: Disposal Considerations

13.1 Waste Treatment Methods Can be disposed of via licensed disposal company

Section 14: Transport Information

- 14.1 UN Number Non-regulated material
- 14.2 UN Proper Shipping Name Non-regulated material
- 14.3 Transport Hazard Class Non-regulated material
- 14.4 Packaging Group Non-regulated material
- **14.5 Environmental Hazards** No
- **14.6** Special Precautions for User Not applicable

Section 15: Regulatory Information

- 15.1 Safety, Health and Environmental Regulations/Legislation Specific for the Substance or Mixture None of the ingredients is listed
- **15.2 Chemical Safety Assessment** A chemical safety assessment has not been carried out

Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate,



BS EN ISO 9001:2015 Certificate No. FM 36099

Registered Office:

Braintree, Essex, CM7 2SE

Registered in England No. 275



Tel: +44 (0) 1376 340 456 Fax: + 44 (0) 1376 340 453 Sales Fax: +44 (0) 1376 340 426 e-mail: <u>sentekuk@btconnect.com</u> website: <u>www.sentek.co.uk</u>

Section 16: Other Information

The information above is based on our present knowledge and is believed to be correct but shall be used only as a guide. It does not constitute a guarantee for any specific product features.

Date of Preparation/Last Revision

22/06/2020



BS EN ISO 9001:2015

Certificate No. FM 36099

Registered Office:

Unit 6-7 Crittall Court, Crittall Drive, Springwood Industrial Estate, Braintree, Essex, CM7 2SE