



SPECIFICATION SHEET

R1 GLASS CALOMEL REFERENCE ELECTRODE

| | |
|------------------------|--------------------------|
| Part Number | 401-75 |
| Temperature Range (C°) | 0 – 50 |
| Cable Length (m) | 1 |
| Connector | 2mm Plug |
| Refillable | Yes |
| Reference Type | Single Junction, Hg/HgCl |
| Junction Material | Frit Ceramic |
| Body Dimensions (mm) | 12 x 120 |
| Body Material | Glass |
| Warranty Off-Shelf | 12 months |

REGISTERED OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

HEAD OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

Tel: +44 (0) 1376 340 456

Fax: +44 (0) 1376 340 453

E-mail: sales@sentek.co.uk

Website: www.sentek.co.uk

Company Registration No. 02753857



FM 36099