



## SPECIFICATION SHEET

### P19/BNC PH ELECTRODE

Part Number	119-75
Body Material	Glass
pH Range	0 - 14
Temperature Range (C°)	0 - 80
Combination	Yes
Refillable	Yes
Cable Length	1 Metre
Connector	BNC
Upper Shaft Length	90mm
Upper Shaft Diameter	12mm
Lower Shaft Length	40mm
Lower Shaft Diameter	6mm
E O Point	0 ± 25mV
Response Time	<30 seconds 95% FSD
Junction	2 x Ceramic Frit
Reference	Ag/AgCl, Liquid Filled, Single Junction
Warranty off the Shelf	12 Months
Cleaning Solution	Cleaning Solution – Order No. 791-15
Storage Solution	3.5M KCl/AgCl Solution – Order No. 551-68
Refill Solution	3.5M KCl/AgCl Solution – Order No. 551-68

REGISTERED OFFICE  
Unit 6-7 Crittall Court, Crittall  
Drive, Springwood Industrial  
Estate, Braintree, Essex,  
CM7 2SE

HEAD OFFICE  
Unit 6-7 Crittall Court, Crittall  
Drive, Springwood Industrial  
Estate, Braintree, Essex,  
CM7 2SE

Tel: +44 (0) 1376 340 456  
Fax: +44 (0) 1376 340 453  
E-mail: sales@sentek.co.uk  
Website: www.sentek.co.uk  
Company Registration No. 02753857

