



SPECIFICATION SHEET

PI10/LAB/BNC PH ELECTRODE

Part Number	140-012
Measuring Range (pH)	0 - 14
Temperature Range (C°)	0 - 130°C (Tip Only)
ATC	N/A
Response Time	Less than 30 seconds, 95% FSD
Cable Length (m)	1
pH Termination	BNC
Combination	Yes
Refillable	No
Reference Type	Double Junction, Gel Filled, Ag/AgCl
Junction Material	Porous Teflon
E O Point (pH 7)	0 ± 25mV
Pressure	10 BAR
Body Dimensions (mm)	12 x 120
Body Material	Glass
Cleaning Solution	Cleaning Solution – Order No. 791-15
Storage Solution	KCl Storage Solution – Order No. 791-65
Warranty Off-Shelf	12 months

REGISTERED OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

HEAD OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

Tel: +44 (0) 1376 340 456
Fax: +44 (0) 1376 340 453
E-mail: sales@sentek.co.uk
Website: www.sentek.co.uk
Company Registration No. 02753857



FM 36099