

SPECIFICATION SHEET

P11/GEL/6X200/PT1000/VP6 PH ELECTRODE

Part Number	111-983
Measuring Range (pH)	0 - 14
Temperature Range (C°)	0 - 80
Response Time	Less than 30 seconds, 95% FSD
Cable Length (m)	N/A
Connector	VP6
Temperature	PT1000
Temperature Connector	VP6
Refillable	No
Combination	Yes
Reference Type	Double Junction, Gel Filled, Ag/AgCl
Junction Material	Ceramic Frit
E O Point (pH 7)	± 25mV
Slope	95 -102%
Body Dimensions (mm)	200 x 6
Body Material	Glass
Thread Size	PG 13.5
Cleaning Solution	Cleaning Solution – Order No. 791-15
Storage Solution	KCl Storage Solution – Order No. 791-65
Warranty Off the Shelf	12 months

REGISTERED OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

HEAD OFFICE

Unit 6-7 Crittall Court, Crittall
Drive, Springwood Industrial
Estate, Braintree, Essex,
CM7 2SE

Tel: +44 (0) 1376 340 456
Fax: +44 (0) 1376 340 453
E-mail: sales@sentek.co.uk
Website: www.sentek.co.uk
Company Registration No. 02753857

